


## AndraBalloon® & Andra2Balloon Stent-delivery Catheter

Optimized for stent dilatation  
Balloon diameter 10 – 32 mm  
Balloon length 20 – 60 mm

# AndraBalloon® & Andra2Balloon

## Stent-delivery Catheter



- Non-compliant balloon made of polyamide (PA)
- Coaxial shaft design for fast inflation and deflation
- Optimized balloon catheter for stent delivery
- 2 X-Ray markers on the inner shaft
- Catheter length: 120 cm (AndraBalloon®); 100 cm (Andra2Balloon)
- Recommended introducer sheath: 7 Fr – 14 Fr
- Guide wire: max. 0.035" ( $\leq$  0.90 mm)

### Order Codes AndraBalloon®

Order Code	Diameter/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
ABL 10/20-1200	10/20	400 (4)	1000 (10)	7
ABL 10/30-1200	10/30			
ABL 10/40-1200	10/40			
ABL 12/20-1200	12/20	400 (4)	1000 (10)	7
ABL 12/30-1200	12/30			
ABL 12/40-1200	12/40			
ABXL 14/30-1200	14/30	400 (4)	1000 (10)	8
ABXL 14/40-1200	14/40			
ABXL 14/60-1200	14/60			
ABXL 16/20-1200	16/20	400 (4)	800 (8)	10
ABXL 16/30-1200	16/30			
ABXL 16/40-1200	16/40			
ABXL 16/60-1200	16/60			
ABXL 18/20-1200	18/20	400 (4)	700 (7)	10
ABXL 18/30-1200	18/30			
ABXL 18/40-1200	18/40			
ABXL 18/60-1200	18/60			

Packaging Unit: 1

Order Code	Diameter/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
ABXL 20/20-1200 ABXL 20/30-1200 ABXL 20/40-1200 ABXL 20/60-1200	20/20 20/30 20/40 20/60	400 (4)	700 (7)	11
ABXL 22/40-1200 ABXL 22/60-1200	22/40 22/60	400 (4)	700 (7)	11
ABXL 24/40-1200 ABXL 24/60-1200	24/40 24/60	400 (4)	700 (7)	12
ABXL 26/40-1200 ABXL 26/60-1200	26/40 26/60	400 (4)	700 (7)	13
ABXL 28/40-1200 ABXL 28/60-1200	28/40 28/60	400 (4)	700 (7)	13
ABXL 30/40-1200 ABXL 30/60-1200	30/40 30/60	400 (4)	600 (6)	13
ABXL 32/40-1200 ABXL 32/60-1200	32/40 32/60	400 (4)	600 (6)	13

Packaging Unit: 1

## Order Codes Andra2Balloon

Order Code	Inner Balloon Ø/Length in mm	Outer Balloon Ø/Length in mm	Nominal Pressure in kPa (bar)	Rated Burst Pressure in kPa (bar)	Introducer Sheath Ø in Fr
A2B 20/40-1000 A2B 20/60-1000	10/40 10/60	20/40 20/60	400 (4)	700 (7)	11
A2B 22/40-1000 A2B 22/60-1000	10/40 10/60	22/40 22/60	400 (4)	700 (7)	12
A2B 24/40-1000 A2B 24/60-1000	14/40 14/60	24/40 24/60	400 (4)	700 (7)	13
A2B 26/40-1000 A2B 26/60-1000	14/40 14/60	26/40 26/60	400 (4)	700 (7)	14
A2B 28/40-1000 A2B 28/60-1000	14/40 14/60	28/40 28/60	400 (4)	700 (7)	14
A2B 30/40-1000 A2B 30/60-1000	14/40 14/60	30/40 30/60	400 (4)	600 (6)	14
A2B 32/40-1000 A2B 32/60-1000	14/40 14/60	32/40 32/60	400 (4)	600 (6)	14

Packaging Unit: 1

Medical Devices – Made in Germany



**Andramed GmbH**

Schiesswieslenstr. 18  
72766 Reutlingen, Germany  
Tel.: +49 (0) 7127 811 57 15  
Fax: +49 (0) 7127 811 57 11  
info@andramed.com  
www.andramed.com

Made in Germany  
Rev. 04/2025-05-06

